Artikel	Articulo
Item	Artikeln
Article	Typ
Artikel	Tovar
Articolo	

4519

4519GAz_1408

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1 | Introduction

1.1 | Theory of therapy

Transcutaneous electrical nerve stimulation (TENS) is a non-invasive technique in which a low-voltage electrical current is delivered through wires from a small power unit to electrodes located on the skin. Electrodes are temporarily attached with paste in various patterns, depending on the specific condition and treatment goals. TENS is often used to treat pain, as an alternative or addition to pain medications. Therapy sessions may last from minutes to hours. The use of electrical stimulus for pain relief was popularized in the 19th century and became widespread in the 1960s and 1970s using battery power.

Transcutaneous electrical nerve stimulation (TENS) was first introduced into current clinical practice following Melzack and Wall's gate control theory of pain in 1965. Davis (1993) and Lewith (1984) explain the gate control theory of pain as follows. An area of the dorsal horn of the spinal cord, known as the substantia gelatinosa, acts as a gate to nociceptive impulses. It receives myelinated nerve fibres, the largest being A fibres, and small non-myelinated nerve fibres (C fibres). If pain impulses pass along fine myelinated fibres and C fibres rather than along A fibres, the gate is opened, and the patient perceives pain. If A fibre transmission of impulses is greater, the gate may be closed.

There is also evidence that the TENS machine enhances the production of the body's own natural pain killing substances: endorphins and enkephalins. Human body produces endorphins and enkephalins, which are opiate-like substances to counter the pain. Low frequency stimulation causes the release of the endorphins and enkephalins.

1.2 | Why consider digital pain relief?

Pain is a warning signal -we need these signals to tell us that something may be wrong with our body. Without it, we may do not know that part of our body might be damaged, thereby damaging them further. However, once we have identified damage, pain serves little purpose. In the case of chronic, regular pain it can significantly interfere with daily activities and the quality of life.

1.3 | How does the Digital Pain Relief (TENS) work??


Digital Pain Relief (TENS) works by passing harmless electrical signals into the body from its pads. This relieves pain in two ways:

- Firstly, it blocks the body's pain signals. These are normally transmitted


from the area of damage through the nerve fibers to the brain, TENS interrupts these pain signals.

- Secondly, TENS stimulates the body's production of endorphins - its own natural painkillers.

2 | Contents and display indicators

- 1 Mainframe
- 2 Batteries
- 3 Instruction manual
- 4 4 Electrodes (Pads)
- 5 2 Wire

On/Off: Press the button to turn on/turn off the unit.


CH1(+): Press the button to increase the CH 1's intensity from 0 to 9.

CH1(-) : Press the button to decrease the CH 1's intensity from 9 to 0.

CH2(+): Press the button to increase the CH 2's intensity from 0 to 9.

CH2 (-): Press the button to decrease the CH 2's intensity from 9 to 0.

Program: Press the button to choose different program in turn (default program as „P1“). Total 16 programs from 1-16.

3 | Intended use

AD-2021 is a TENS device which is effective in relieving pain.

Treatment effect:

- Relieving lower-back pain
- Stimulate muscles
- Promoting blood circulation
- Eliminate tiredness

4 | Contraindication

The device is safe for most all people, with the following exceptions:

- People with acute disease
- Cancer patients
- People with infectious skin wounds
- Women who are menstrating or expectant mothers
- People with heart disease
- People with high fever
- People with abnormal blood pressure
- People who have no feeling about their skins or people with abnormal skins
- People receiving physiotherapy or with unusual sensations of their body other than noted above.

5 | Product description

The AD-2021 is a battery operated pulse generator that sends electrical impulses through electrodes to the body and reaches the underlying nerves and muscle groups. The device is provided with two controllable output channels, each independent of each other. An electrode pair can be connected to each output channel.

The electronics of the AD-2021 Digital create electrical impulses whose Intensity, Pulse Width, Pulse Rate may be altered according to the program. Press buttons are very easy to use and the panel cover prevents changes in the setting.

The AD-2021 corresponds to the below standards:

- „IEC 60601-1: 2005 + CORR. 1 (2006) + CORR. 2 (2007)/EN 60601-1:2006/AC:2010 (Medical electrical equipment -- Part 1: General requirements for basic safety and essential performance)“
- IEC60601-1-2:2007/EN 60601-1-2:2007 /AC:2010 (Medical electrical

equipment -- Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests)

6 | Specifications

Product name:	TENS Device - 2 canals
Model:	AD-2021
Number of treatment program:	16
Range of pulse frequency:	2 ~ 50 Hz
Range of output voltage:	30 V \pm 20%(500 Ohm load)
Classification:	Internally powered, Type BF applied part, IPX0, No AP or APG, Continuous operation
Machine size:	Approx. 120mmx62mmx19mm
Weight:	Approx. 87g (excl. battery)
Power source: batteries:	2 x 1.5 V SIZE AAA
Environmental temperature	
- for operation:	-10°C - 40°C
- for storage and transport:	-20°C - 55°C
Environmental humidity	
- for operation:	<80%
- for storage and transport:	<95%
Environmental pressure:	80 KPa - 105 KPa
Battery life:	Approx. 2 months with alkaline batteries and 30 min. usage per day.

7 | Therapy mode

NO.	PROGRAM	APPLICATION	TIME
1	Light muscle stimulation	Initial support of chronically injured muscle	45 - 60 Min
2	Aerobic endurance	Training for muscle groups targeted by aerobic sports	10 - 15 Min
3	Massage	Relaxing, restful warming, stimulates capillaries	45 - 60 Min
4	Muscle stimulation	Stomach and Thighs	45 - 60 Min
5	Muscle stimulation	Stomach, Legs, Arms	30 - 45 Min
6	Profusion/Lymphatic drainage	Improves localized circulation and profusion	45 - 60 Min
7	TENS Contractions	Limbering effect / muscle pain	45 - 60 Min
8	TENS Contractions	Accute muscle pain / cramping	45 - 60 Min
9	Muscle stimulation	Standard massage	30 - 45 Min
10	Muscle stimulation	Fitness massage	10 - 15 Min

11	Muscle stimulation	Enduring massage	10 - 15 Min
12	Muscle stimulation	Skin massage	45 - 60 Min
13	Muscle stimulation	Contraction	10 - 15 Min
14	TENS Relax 1	Relieves fatigue at the end of the day	45 - 60 Min
15	TENS Relax 2	Signs of fatigue/ pain	45 - 60 Min
16	Post accident Rehabilitation	Muscle stimulation to rehabilitate	45 - 60 Min

8 | Notice

- If you feel uncomfortable or abnormal skin sensations, please stop using. Consult with doctors and follow doctor's advice.
- Before you are to shift the electrode pad to a different position in using the apparatus, you must turn off the power first.
- Please don't use it in the bathroom or other place with high humidity.
- Please don't use it when driving or sleeping.
- Do not make any sharp kinks in the connecting leads or electrodes.
- Please don't use it for purposes other than outlined in this manual.
- In the process of stimulating and therapy, please do not allow metal parts of belts, wristwatches or necklaces touch the leaf electrode pads.
- Do not use the device if you are connected to, or in the vicinity of, high-frequency surgical or industrial equipment. This may cause burn injuries on the skin under the electrodes, as well as problems with the stimulator.
- Observe caution when using the device in the immediate vicinity of cellular phones that are switched on.
- Do not use the device in the direct vicinity of short-wave or microwave equipment, since this may affect the output power of the stimulator.
- This TENS Device is designed for adults and should never be used on infants or young children. Consult your physician or other health care professionals before use on older children.
- The device might not meet its performance specifications or cause safety hazard if stored or used outside the specified temperature and humidity ranges in this manual.
- Please do not place pads near heart, genitals or diseased skin.
- Users with implanted electronic equipment, such as pacemakers and intracardiac defibrillators who have not gotten doctor's permission must not use the device. Pregnant women should not use the device during the first trimester, and should always consult a doctor, midwife or physiotherapist prior to use.
- Simultaneous connection of a PATIENT to a h.f. surgical EQUIPMENT may result in burns at the site of the STIMULATOR electrodes and possible

damage to the STIMULATOR.

- Operation in close proximity (e.g. 1 m) to a shortwave or microwave therapy EQUIPMENT may produce instability in the STIMULATOR output.
- Application of electrodes near the thorax may increase the risk of cardiac fibrillation.
- Please do not drop, repair, or rebuild it privately.
- Please do not use electrode pads other than those supplied by the manufacturer, otherwise it may bring biocompatible hazard and may result in errors.
- Please do not share the electrode pads with infective persons to avoid cross-infection.

9 | Setup and operating procedures

1. Battery loading

- a Open battery cover at the back of the device.
- b Load two „AAA“ size batteries. Please pay attention to polarity.
- c Close the battery cover.

- Rechargeable batteries are not suitable for this device.
- Remove the batteries if the monitor will not be used for a month or more to avoid relevant damage of battery leakage.
- Please do not mix use new and old batteries or different type of batteries.
- Please do not throw battery into fire.

Disposal of appliances and batteries

The device and the batteries, must be disposed of according to local regulations at the end of their usage.

2. Attachment of electrode lead wires


- a Insert the plug into the electrode socket on the main apparatus.
- b Connect the snap-fasteners at the other end of the wire to the electrode pads.
- c Place the electrode on your body as directed by your physician.


- Clean the wires by wiping with a damp cloth.
- Coating them lightly with talcum powder will reduce tangling and prolong life.

3. Operating instructions


- a Turn on the „On/Off „button.

Increase
„Intensity“ for
CH1 (CH1)

Decrease
„Intensity“ for
CH1 (CH1)

Increase
„Intensity“ for
CH2 (CH2)

Decrease
„Intensity“ for
CH2 (CH2)

- b When one of two channels has not been connected appropriately, it will not be displayed on the LCD.
- c Press the „Program“ button to choose different program in turn (default program as „P1“), total 16 programs from 1-16.
- d Press the „Intensity“ button of each output channel to choose different intensity. (default intensity as „0“), total 9 intensities from 0-9.
- e Finishing the Operation










Pull out the plug from the main apparatus. Clean the electrodes.

You should hold the plug when pulling it out. Please do not pull the wire.

Clean the electrodes:

- You should hold the plug when pulling it out. Please do not pull the wire.
- Please use clear water to wash or use wet cloth (instead of facial tissue) to gently wipe up the electrodes when cleaning them. Do not use brush or fingernail to clean surface of the electrodes to avoid scratching them.

4. Usage diagrams

Upper back**Lower back****Shoulder****Chest****Belly****Arm****Adductor and gluteus****Legs****Feet****Joint**

10 | Troubleshooting

PROBLEM	POSSIBLE CAUSE	SOLUTION
You have no feeling of stimulus.	Are the batteries exhausted?	Replace the batteries.
	Are the batteries properly loaded?	Correctly load the batteries.
	Is the wire properly connected?	Firmly connect the wire.
	Have you removed the transparent protection film over the electrode pad?	Remove the transparent protection film over the electrode pad.
Stimulus is weak.	Do the electrode pads firmly stick to the skin?	Firmly stick the electrode pad to the skin.
	Are the electrode pads overlapped?	Separate the electrode pad and stick them to the skin gain.
	Are the electrode pads dirty?	Please clean the electrode pad.
	Is intensity too weak?	Turn the intensity regulation dial to regulate it.
	Are the electrode pads properly positioned?	Change the position of the electrode pad.
The skin becomes red.	Is the therapeutic time too long?	Control it within 10-15 minutes a time.
	Do the electrode pads firmly stick to the skin?	Please firmly stick the electrode pad to the skin.
	Are the electrode pads too dry?	Please gently wipe them up with wet cloth and then use them again.

The skin becomes red.	Are the electrode pads dirty?	Please clean the electrode pad.
	Are the surface of the electrode pads scratched?	Please replace them with new electrode pad.
Power sources cut off in the therapeutic process.	Have the electrode pads come off the skin?	Turn off the power and stick the electrode pad firmly to the skin.
	Are the wires disconnected?	Turn off the power and properly connect the wire.
	Have the batteries been exhausted?	Please replace them with new ones.

11 | Maintenance

- Do not drop this device subject it to strong impact.
- Avoid high temperature and direct sunlight. Do not immerse the device in water.
- If this device is stored near freezing, allow it to acclimate to room temperature before use.
- Do not attempt to disassemble this device.
- If you do not use the device for a long time, please remove the batteries.
- If the device becomes dirty, please clean it with a soft dry cloth. Do not use any abrasive or volatile cleaners.
- No component can be maintained by user in the device. The circuit diagrams, component part lists, descriptions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated as repairable can be supplied by the manufacturer.
- In order to ensure that the electrodes adhere as long as possible, they should be cleaned carefully with a damp, lint-free cloth.
- After use, stick the electrodes onto the backing film.
- Electrode Pads are designed to be replaceable, so when the electrodes do not adhere firmly, you must replace electrode-pads.

12 | Explanation of symbols on unit

Symbol for „CONSULT Operation Guide“.

Symbol for „WARNING!“

Symbol for „TYPE BF APPLIED PARTS“.

Environment Protection — Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local Authority or retailer for recycling advice.

Symbol for „MANUFACTURER“

Symbol for „DATE OF MANUFACTURE“

Symbol for „KEEP DRY“

SN

Symbol for „SERIAL NUMBER“

CE 0483

Complies with MDD93/42/EEC requirements

13 | Warranty information

The guarantee period is 24 months from the date of purchase and covers all manufacturers' errors in material and quality. This guarantee only applies if the instructions for use have been followed and becomes ineffective if force has been applied to the appliance or it has been used in any but the correct and proper manner, or if it has been serviced by any unauthorised person.

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Hydas GmbH
c/o atrikom fulfillment GmbH • Haagweg 12
Hydas GmbH & Co.KG
Andon Health Co., Ltd.
Hydas Inc.
Hummelstown, PA. 17036 U.S.A.